



T-REX Series

OPHARDT
hygiene



New

**HAND HYGIENE
MADE EASY!**

RX 5-T

Desktop touchless dispenser

A mobile, flexible dispensing solution providing hand hygiene opportunities to increase compliance standards.

- 500 ml-Euro bottle compatible 
- Ready for disposable hygiene packaging
- Optimized footprint where space is limited
- For use with disinfectant and soaps 
- Free standing or wall mounted
- 4 C cell battery operated
- 1.5 ml dosage



Hygiene stand
ordered separately

